

4/8/2022 - 4/12/2022

FRIDAY, 4/8/2022

Day 1 Objective

- Our goal today is to amplify the upstream and downstream regions of the ACIAD2049 ADP1 gene
- We also identified a gene within ADP1 to delete and use the genome to design a primer for that gene.

PCR Reactions

Thermocycler Settings

- 50 μ L
- Denaturing- 98
- Annealing- 58
- Elongation- 72

- 10 μ L of Phusion HF buffer
- 1 μ L of mM dNTPs
- 2.5 μ L Forward Primer
- 2.5 μ L Rverse Primer
- .5 μ L Phusion DNA polymerase
- 1.5 μ L DMSO
- 31 nuclease free water
- 1 μ L template DNA

Primer Design

Golden Gate

gc 40-60%

melting temp 55-60

forward and reverse withing 5 degrees of each other

Forward part creation primer 5' extension

SATURDAY, 4/9/2022

Day 2 objective-

- 1.) First we are running gel electrophoresis on each of the 6 PCR products to determine if DNA was successfully made and if we made the correct dna sequence
- 2.) Finally, we're eluting DNA from the PCR products with the PCR cleanup kit and then determining the [DNA] with the Quibit

Agarose Gel Electrophoresis

- Dissolve agarose in TE buffer
- Add SYBR Safe dye
- Set Up casting tray with 15 well comb
- Pour molten agar into tray, allow to solidify for 25-30 minutes
- Transfer gel into electrophoresis unit
- Load ladder into one lane of the gel
- Mix 1 μ L of dye with 5 μ L of PCR product. Load all 6 products into a lane
- Make sure to label each lane
- Start electrophoresis running at 130 volts for 20 minutes\
- Visualize gel, upload to benchling

PCR Cleanup-

- Follow instructions in DNA clean and Concentration kit.
- For elution step, add 20 μ L of buffer, incubate at room temp for 1 min befor centrifuging

-Quantify DNA of final elution with qubit (see week 1 procedure)

DNA Concentrations-

TUESDAY, 4/12/2022

Colony Counts

DNA Control 1 10^{-2} LB+AR+KAN

- Positive Control- 1024
- Negative Control- 0
- Transformation Efficiency = $(1024+0)/2 * 10^6 = 5.12E8$

DNA Control 1 10^{-5} LB+AR

- Positive Control- 556
- Negative Control- 540
- Transformation Efficiency- $(556+540)/2 * 10^6 = 5.48E8$

DNA Control 2 10^{-2} LB+AR+KAN

- Positive Control- 608
- Negative Control 0
- Transformation Efficiency- $(608)/2 * 10^6 = 3.04e8$

DNA Control 2 10^{-5} LB+AR

- Positive Control- 444
- Negative Control- 512
- Transformation Efficiency- $(512+444)/2 * 10^6 = 4.78E8$

2.2e9 for tdkkan