

13 uL dH20Follow Up-

o Run a PCR of this assembley to amplify the entire sequence

ABR Detection Threshold

- Grew up a culture of 013 strain for detection threshold
 - Follow up-
 - Perform a transformation with the decreasing DNA concentrations of NPTII
 - Perform a transformation with the decreasing DNA concentrations of Tem-1


Co-Culture Experiment

- Grew up a culture of 019 for TDK Transformation into 2049
 - Follow Up-
 - Add 35 uL of the overnight culture to 500 uL LB + 2049 tdk/kan sequence
- Ran a PCR of the 2049 Tdk/kan gene
 - o 001 AND 004
 - o Run at 59° with an extension time of 2:15

TUESDAY, 8/23/2022

Co-Culture Experiment-

- Ran the gel for the 2049 PCR I ran yesterday
- Gel Order: Ladder Control tdk tdk GGA



P. Destructans Detector Test-

- Run a PCR pf yesterday's Golden Gate Assembley of the two PCR's
 - $\circ\quad$ This PCR will ligate the two PD sequences to make the entire target sequence
- Primers = 063 + 066
 - o Run at 64°
 - o Run at 66°
 - ✓ 124 uL of water
 - 40 uL of HF buffer
 - ✓ 10 uL of 063
 - ✓ 10 uL of 066
 - 6 uL DMSO
 - 4 uL dntp
 - 2 uL polymerase

2049 Up + PD 3a (1200+500) 30 sec /kb 001+064 2049 Repressor + PD 3 h + 2049 down (1000+500+1200)
004 + 061

Print Settings

Paper Size

Letter (8.5" x 11")

Layout

Fit To Screen
Portrait Landscape

Customization
Show link URLs

Print