

Technic : Transformation

Goal: Introduce recombinant plasmid DNA (ligation product) into competent *E. coli* XL10 Gold by heat shock, followed by recovery and plating on selective medium.

Procédure:

1. Thaw competent cells. Retrieve XL10 Gold competent cells from the -80°C freezer and thaw on ice.
2. DNA addition. Add $10\text{ }\mu\text{L}$ of ligation product to $100\text{ }\mu\text{L}$ of thawed competent cells in a sterile tube. Mix gently by flicking; avoid pipetting up and down.
3. Incubation on ice. Keep the mixture on ice for 30 min at 4°C .
4. Heat shock. Transfer the tube to a 42°C water bath or heat block for 30 seconds.
5. Recovery on ice. Immediately place the tube back on ice for 2 min.
6. Add $150\text{ }\mu\text{L}$ of LB
7. Outgrowth. Incubate the tube at 37°C with shaking (200–250 rpm) for 1 hour to allow expression of antibiotic resistance genes.
8. Plate the resuspension onto a prewarmed LB agar plate containing ampicillin. Spread evenly across the surface ("like a crêpe") to ensure well-distributed colonies.
9. Drying. Allow the plate surface to dry briefly at room temperature until no liquid remains visible.
10. Incubation. Place the plate inverted at 37°C overnight to allow colony growth.

Reagents and Solutions:

- XL10 Gold chemically competent cells (stored at -80°C)
- Ligation product (DNA of interest)
- LB agar plates supplemented with antibiotic (ampicillin, final concentration $100\text{ }\mu\text{g/mL}$)

Equipment:

- Sterile microcentrifuge tubes
- Ice bucket
- Water bath or heat block set to 42°C
- Incubator/shaker set to 37°C
- Sterile spreaders or glass beads
- Pipettes and sterile tips
- Petri dish rack

Estimated time: ~ 3 hours

Controls and Troubleshooting:

- Verification of bacterial vitality on antibiotic-free agar, if no bacterial growth in selective media, perhaps antibiotic resistance hasn't had time to develop. In this case, incubate longer at step 7.
- Verification of plasmid antibiotic resistance by spreading untransformed bacteria on selective agar