

# Lab Notebook

## Protocols from papers:

- Protocol for DNA extraction from PCR diagnostics paper (Toda et al, 2014): Adult females were homogenized individually in 10  $\mu$ l of lysis buffer including 9  $\mu$ l of STE buffer [100 mM NaCl, 1 mM EDTA (pH 8.0), 10 mM Tris-HCl (pH 8.0)], and 1  $\mu$ l of proteinase K (10 mg/ml). The homogenate was incubated for 15 min at 65 °C and then heated for 5 min at 95 °C to inactivate the proteinase K . PCR was per-formed in 20  $\mu$ l reaction volume, containing 20 pmol of each primer, 10 mM Tris-HCl (pH 8.3), 50 mM KCl, 2.5 mM MgCl<sub>2</sub>, 0.25 mM of each dNTP, and 0.5 units of Taq-polymerase (AmpliTaq Gold, Applied Biosystems). The PCR program in the DNA thermal cycler (DNA thermal cycler 2720, PE Applied Biosystems) was 94°C for 9 min as a preheat step, followed by 35 cycles at 94°C for 30 s, 52°C for 30 s, and 72°C for 1 min, and then 72 °C for 7 min as a final extension.

- Protocol for DNA extraction from LAMP diagnostics paper (Toda et al, 2015): Place one individual *Scirtothrips dorsalis* insect into a 0.2 ml PCR tube (using a fine-tipped brush). Add 10  $\mu$ l of 100 mM Tris-HCl (pH 8.0) to the tube. Grind the insect thoroughly in this solution (using a fine-tipped brush or a heat-modified pipette tip as mentioned previously). Heat the homogenate at 99.9°C for 10 minutes. The resulting supernatant/lysate containing the DNA is then used directly as the template for the LAMP reaction (0.5  $\mu$ l of this extract was used per 20  $\mu$ l LAMP reaction).

## **2025/05/04. Field sampling**

- Visit to the Mango farm.
- Collected 6 eppendorf tubes of thrips (approximately >10 thrips in each tube). They are stored in 99% ethanol.

## **2025/06/11 DNA extraction of collected samples.**

Quick Protocol: DNeasy Blood & Tissue Kit (for insects)

Materials Needed:

- DNeasy Blood & Tissue Kit (Qiagen)
- 1.5 mL microcentrifuge tubes
- Ethanol (96–100%)
- Proteinase K
- Buffer ATL, AL, AW1, AW2, AE

---

### 1. Sample Prep

- 5 samples collected from farms were used (4 samples collected by Baha and 1 sample collected by Wahei)
- Samples were centrifuged to spin down the insects and ethanol was discarded
- Incubated in 56 °C until excess ethanol dried out
- Add 180 µL Buffer ATL.
- Add 20 µL Proteinase K.
- Vortex briefly or flick.
- Incubate at 56 °C for 1–3 hours (or overnight) until fully lysed.

### 2. Lysis

- Add 200 µL Buffer AL, mix.
- Add 200 µL ethanol (96–100%), mix well.

### 3. Binding

- Transfer the mixture to a DNeasy spin column in a collection tube.
- Centrifuge 1 min at  $\geq 6000 \times g$  (8000 rpm).
- Discard flow-through.

### 4. Wash

- Add 500 µL Buffer AW1, spin 1 min, discard flow-through.
- x

### 5. Elution

- Place column in a clean 1.5 mL tube.
- Add 100 µL Buffer AE (or nuclease-free water).

- Incubate 1 min, then centrifuge 1 min at  $\geq 6000 \times g$ .

Store:

- Store eluted DNA at  $-20^\circ\text{C}$ .

**2025/06/12 PCR tests of extracted DNA samples.**

PCR reaction mix prep:

- 12.5 ul of Phusion Plus Master Mix
- 1.25 ul of 10uM of each ThripsITSF
- 1.25 ul of 10uM of each ThripsITSR
- 1.25 ul of 10uM of each SdITSF3
- 1.25 ul of 10uM of each SdITSF7
- 3 ul of extracted DNA
- 4.5 ul of H<sub>2</sub>O

PCR conditions:

Initial denaturation	98°C	30 sec	
Denaturation	98°C	10 sec	} 25 cycles
Annealing	60°C	10 sec	
Extension	72°C	30 sec	
Final extension	72°C	5 min	

Run on 1% agarose gel (0.5 g of agarose and 50 ml of TAE buffer) with 180 V for 5 min

The result was unclear.

## 26/06/2025 PCR test extracted DNA samples (optimize)

This time, we used PCR protocol from paper:

PCR reaction mix prep:

- 12.5 ul of Phusion Plus Master Mix
- 1 ul of 20uM of each ThripsITSF (20 pmol)
- 1 ul of 20uM of each ThripsITSR (20 pmol)
- 1 ul of 20uM of each SdITSF3 (20 pmol)
- 1 ul of 20uM of each SdITSF7 (20 pmol)
- 5 ul of extracted DNA
- 3.5 of H2O


PCR conditions:

Initial denaturation	94°C	9 min	
Denaturation	94°C	30 sec	} 35 cycles
Annealing	52°C	30 sec	
Extension	72°C	1 min	
Final extension	72°C	7 min	

Mix PCR samples with loading dye (5:1) and use the appropriate ladder (with bands on 200-500 bp). Run on gel 2% agarose gel (1 g of agarose dissolved in 50 ml of TAE) with 130 V for 15 min.

Clear bands near 200 bp (characteristic to C strain), one sample has band on 300 bp (characteristic to 300 bp).

Samples were loaded in this order: diluted ladder, miss, 7ng/ul extracted dna (trap water), 3.3 ng/ul, 2.6 ng/ul, 0.8 ng/ul, 0 ng/ul (Figure 1).



### 13/07/2025- Sticky sheet DNA extraction optimization

**Aim:** To extract DNA from the water sample where the sticky sheet was soaked

#### A. Material

250713

Sticky Sheet in the RO water (Kept in my room)

Sticky Sheet in the 80% ethanol (Kept in 2 tube in -20 freezer)

Leaf and thrips in the 80% ethanol (Kept in 4 tube in -20 freezer)

Thrips in the 80% ethanol (Kept in 4 tube in -20 freezer)

3 water bottles

#### B. Procedure

1. Get sticky sheet from farmer (250713)
2. Soak it in the bottle filled with 500ml water (250713)
3. Take 5 ml water from the bottle and put them in the tube everyday (from day0-day7)  
Day0: 250713
4. Freeze the tube
5. Extract DNA from these tubes in same procedure in previous experiment
6. Check PCR amplification

## **2025/07/29 DNA extraction from Stick sheet sample (SS1)**

Wahei and Youjung

### **Aim**

To extract DNA from the samples for 7 days where the sticky sheet was soaked by the protocol (**2025/06/11**)

### **Procedure**

#### **Water sample: 200 µl**

See protocol 2025/06/11 DNeasy Blood & Tissue Kit (for insects)

Conduct from “Add 180 µL Buffer ATL.”

### **Result**

Get extraction from 7 samples

### **Next**

They will be checked by PCR following methods on **26/06/2025** in 08/01

### **Rest of work**

#### **1. Extraction dna from water sample**

Sticky Sheet (SS) (n=3)

7 sample for 7 days (or 14 days) each replicate

Total reaction 63 reactions for PCR

SS1

- PCR on extracted DNA on 250729
- PCR on water
- PCR on heating water

SS2 and SS3

- Extract DNA and PCR
- PCR on water
- PCR on heating water

Thrips water (T) (n=2)

7 sample for 7 days each replicate (Not yet)

T1 and T2

- PCR on water and heating water
- Extract DNA and PCR

After above, Conduct qPCR on positive samples

## 2025/08/01 PCR

Sticky Sheet (SS) (n=3)

8 sample for 8 days (from day0)

Total reaction 24 reactions for PCR

SS1

- PCR on extracted DNA on 250729
- PCR on water
- PCR on heating water

DNA Concentration

- Water sample

Day0 2.1ng/ul Day1 2.7ng/ul Day2 3.0ng/ul Day3 3.6ng/ul Day4 4.3ng/ul Day5  
4.8ng/ul

Day6 5.3ng/ul Day7 5.8ng/ul

- Extracted sample

Day0 1.8ng/ul Day1 1.3ng/ul Day2 3.3ng/ul Day3 0.0ng/ul Day4 0.1ng/ul Day5  
1.2ng/ul

Day6 1.1ng/ul Day7 0.4ng/ul

PCR reaction mix prep (250625)


- 12.5 ul of Phusion Plus Master Mix
- 1 ul of 20uM of each ThripsITSF (20 pmol)
- 1 ul of 20uM of each ThripsITSR (20 pmol)
- 1 ul of 20uM of each SdITSF3 (20 pmol)
- 1 ul of 20uM of each SdITSF7 (20 pmol)
- 5 ul of extracted DNA
- 3.5 of H2O

PCR conditions:

Initial denaturation	94°C	9 min	
Denaturation	94°C	30 sec	} 35 cycles
Annealing	52°C	30 sec	
Extension	72°C	1 min	
Final extension	72°C	7 min	



Mix PCR samples with loading dye (5:1) and use the appropriate ladder (with bands on 200-500 bp). Run on gel 2% agarose gel (2 g of agarose dissolved in 100 ml of TAE) with 130V for 15 min.

Samples were loaded in this order: ladder, positive control, extracted DNA (day1-7, day0), heated samples (day1-7, day0), water samples(day1-7, day0). No bands.



**PCR of water samples.**

**E1 - extracted samples from water with thrips replicate #1**, E2 - extracted samples from water with thrips replicate #2. H1 - heated samples from water with thrips replicate #1, H2 - heated samples from water with thrips replicate #2, W1 - crude water samples with thrips replicate #1, W2 - crude water samples with thrips replicate #2



## 2025/08/06 eDNA Extraction from Ethanol-Preserved Environmental Samples

**Method:** Qiagen DNeasy Blood & Tissue Kit

**Sample types:** Ethanol-preserved leaf wash & soil samples

**Batch size:** 4 samples

**Total time:** ~1 hr 45 min

---

heat block/incubator (56°C)

---

### Pre-Preparation (Ethanol Removal & Pellet Concentration)


1. **Centrifuge 50 mL Falcon tubes at 4,000 × g for 20 minutes at room temp** to pellet debris.
2. **Carefully decant ethanol** without disturbing pellet.
3. **Air-dry pellet** with lids open in clean area (or 37°C incubator) for **5–10 minutes**.

### DNA Extraction Protocol


5. **Add lysis solution to each pellet:**
  - 180 µL Buffer ATL
  - 20 µL Proteinase K  
(*double volumes if pellet is large*)
6. **Vortex briefly** or pipette up/down to mix.
7. **Incubate at 56°C for 1 hour.**
8. After incubation, add:
  - 200 µL Buffer AL (mix well)
  - 200 µL 96–100% ethanol (mix again)
9. **Transfer lysate (~600 µL total)** to DNeasy spin column in 2 × 300 µL steps.
  - Centrifuge at **6,000 × g (8,000 rpm) for 1 min**, discard flow-through.
10. **Wash 1:**
  - Add 500 µL Buffer AW1
  - Centrifuge 6,000 × g for 1 min
  - Discard flow-through
11. **Wash 2:**
  - Add 500 µL Buffer AW2
  - Centrifuge at **20,000 × g (14,000 rpm) for 3 min**
  - Discard flow-through
12. **Dry spin:**
  - (Optional) Repeat spin for 1 min to ensure all buffer is removed
13. **Elute DNA:**
  - Add 100 µL Buffer AE directly to membrane
  - Incubate 1–2 min at room temp

- Centrifuge at  $6,000 \times g$  for 1 min into clean 1.5 mL tube
- Optional: re-elute with same or new 100  $\mu$ L

From Onna village farm: 3 leaf samples and one soil sample preserved in ethanol  
 From Nanjo farm: 3 leaf samples preserved in ethanol and 3 Q-tip samples (rubbed on different leaves) preserved in AE buffer  
 From Trafru farm: 3 leaf samples preserved in ethanol and 3 Q-tip samples (rubbed on different leaves) preserved in AE buffer  
 Results: 18S: Very good bands from all samples => we have 17 samples to send for sequencing, including the PCR blank



COI: We obtained bands from the samples that were preserved in AE buffer only. So basically, 3 samples from Nanjo farm (good bands), and 3 samples from Trafru (weak bands) => we have 7 samples to send for sequencing including the PCR blank  
 \*Concern: no samples from Onna village farm were preserved in AE buffer so we don't have any samples that amplified with COI from this farm




**2025/09/23 LAMP reaction using thrips DNA.**


Prepared a reaction mix as follows:

<b>Component</b>	<b>25 <math>\mu</math>l Reaction</b>	<b>Final Conc</b>
10X ThermoPol Buffer	2.5 $\mu$ l	1X (contains 2 mM MgSO <sub>4</sub> )
MgSO <sub>4</sub> (100 mM)	1.5 $\mu$ l	6 mM (8 mM total)
dNTP Mix (10 mM)	3.5 $\mu$ l	1.4 mM each
FIP/BIP Primers (25X)	1 $\mu$ l	1.6 $\mu$ M
F3/B3 Primers (25X)	1 $\mu$ l	0.2 $\mu$ M
LoopF/B Primers (25X)	1 $\mu$ l	0.4 $\mu$ M
Bst DNA Polymerase, Large Fragment (8,000 U/ml)	1 $\mu$ l	320 U/ml
DNA Sample	1 $\mu$ l (2-3 ng)	> 10 copies or more
Either phenol red or hydroxy naphthol blue	1-5 $\mu$ l	
Nuclease-free Water	to 25 $\mu$ l	
Total Reaction Volume		25 $\mu$ l

Incubate at 63 C for 60 min. Here is the result after the incubation.



Run samples on a gel with SybrSafe. 1st lane - LAMP reaction without a dye, 2st - Negative control without a dye, 3rd - LAMP reaction with hydroxy naphthol blue, 4th - Negative control with hydroxy naphthol blue.



For phenol red no color change was observed.

Also load the 10 ul of reaction mix with DNA onto Whatman Grade 1 paper squares of ~5mm length and width, however paper quickly absorbs the reaction so color is very bleak.

## PCR with water samples collected by blue light trap

Sample: 2 water samples from blue light trap, 2 insect samples, positive control, negative control

Extract DNA from 2 insect samples by DNeasy Blood & Tissue Kit according to the protocol (25/06/11)


PCR reaction mix prep (same as 25/06/25)

- 12.5 ul of Phusion Plus Master Mix
- 1 ul of 20uM of each ThripsITSF (20 pmol)
- 1 ul of 20uM of each ThripsITSR (20 pmol)
- 1 ul of 20uM of each SdITSF3 (20 pmol)
- 1 ul of 20uM of each SdITSF7 (20 pmol)
- 5 ul of each sample
- 3.5 of H2O

PCR conditions (same as 25/06/25):

Initial denaturation	94°C	9 min	
Denaturation	94°C	30 sec	} 35 cycles
Annealing	52°C	30 sec	
Extension	72°C	1 min	
Final extension	72°C	7 min	

Gel image:



First is positive control. Second lane and third lane is water sample #1 and #2. Fourth and fifth lane are extracted DNA from insects collected by blue light LED trap #1 and #2.

Result: different bands for trap water samples and insect DNA extraction system from positive control indicating 200 and 500 bp bands.