
iGEM - Lysogen Efficiency Testing

Introduction

This protocol is used to determine the degree of lysogen of a temperate phage. It was adapted from the SEA-PHAGES discovery program from the uOttawa TMM teaching laboratories.

Materials

- › Phage lysate
- › *Staph epi 71* Culture
- › TGY plates
- › Phage buffer

Procedure

Day 1

1. Prepare bench for aseptic work and assemble supplies.
 2. Dilute the phage lysate to 1×10^9 pfu/mL, in a final volume of 500 μ L of phage buffer.
 3. Seed plates with phage:
 - a. Obtain 5 agar plates and label them with the name of the phage, the name of the host bacterium (*S. epi 71*), and one of the following dilutions 10^{-2} , 10^{-3} , 10^{-4} , 10^{-5} , and 10^{-6} .
 - b. To each of these 5 plates, add 100 μ L of your phage sample that has been diluted to 1×10^9 pfu/mL.
 - c. Spread the lysate evenly across the plate using sterile glass beads.
 4. Prepare a dilution series of the host bacteria:
 - a. Arrange 5 microcentrifuge tubes in a rack and label them 10^{-1} , 10^{-2} , ..., 10^{-6} .
 - b. Fill each microcentrifuge tube with 450 μ L of phage buffer.
 - c. Add 50 μ L of your undiluted host bacteria culture to the tube labelled " 10^{-1} " and mix by pipetting up and down, and vortexing.
 - d. Using a new clean pipette tip, transfer 50 μ L of the " 10^{-1} " sample to the " 10^{-2} " tube and mix well.
 - e. Continue each successive dilution until you get to your last tube.
 5. Spread host bacteria on control plates:
 - a. Obtain 5 new agar plates and label them with the name of the host bacterium (*S. epi 71*), and one of the following dilutions 10^{-2} , 10^{-3} , 10^{-4} , 10^{-5} , and 10^{-6} . Note: these are your control plates.
 - b. To each of the 5 new agar plates, add 100 μ L of host bacteria at corresponding dilution. Ex; add 100 μ L of host bacteria at the 10^{-2} dilution onto the plate labelled with that dilution.
-

-
- c. Spread the culture evenly across the plate, using sterile beads.
6. Spread host bacteria on phage-seeded plates:
 - a. Obtain the 5 agar plates that have been already seeded with phage.
 - b. To each of these 5 plates, add 100 μ L of host bacteria at corresponding dilution. Ex; add 100 μ L of host bacteria at the 10^{-2} dilution onto the plate labelled with that dilution.
 - c. Spread the culture evenly across the plate, using sterile beads.
 7. Incubate the plates at 34°C until colonies form. For *S. epi*, should take 1-2 nights.

Day 2-3

8. Check plates each day for colony formation.

Colonies that form on the phage-seeded plates are mostly lysogens, though a small fraction of host bacteria may contain a mutation that renders them insensitive to infection by your phage and can also form colonies on the phage-seeded plates.
 9. Once colonies are clearly visible, photograph all plates.
 10. Calculate the rate of lysogeny, also known as the Efficiency of Lysogeny (EOL).
$$\text{EOL} = [(\text{CFU on a phage-seeded plate}) / (\text{CFU on the corresponding control plate})] \times 100$$

Note: CFU = colony forming units. The best plates to use for counting colonies contain 20-500 colonies. The efficiency of lysogeny for a temperate phage can range from 0.1% to 100%.
-